




*Wellness
in Dentistry*



XWELLbed




XWELLbed

The beauty of simplicity

At Swident, we believe dental care should feel effortless for both patients and professionals. That's why we create equipment that enhances comfort, inspires confidence, and makes every treatment more efficient.

The **XWELLbed** embodies this philosophy. By stripping away the unnecessary, we've crafted a design that is pure, streamlined, and purposeful - maximizing efficiency while minimizing risks.

Every movement feels natural and fluid, supported by carefully refined features chosen to empower accurate diagnoses and effective treatments. With **XWELLbed**, simplicity becomes the highest form of sophistication.




Ergo Concept for Dentistry

Design that protects your passion

Dentistry is a lifelong profession, and ergonomics is the key to sustaining it. Our products are built to support a natural posture and relaxed movements, so dentists can work comfortably and confidently for years to come.

Efficiency doesn't just benefit patients with faster, more precise treatments - it also protects the practitioner's body, helping careers last longer. With **XWELLbed**, every detail is designed to keep you in an optimal position, preserving both precision and well-being.



Natural, Relaxed Posture

Swident equipment begins with people. First, we create space that allows practitioners to move freely and naturally. Only then do we shape the tools - trays, angles, and configurations - to serve that freedom. The result is a workspace where accuracy and speed thrive, while fatigue fades away.




Four-Hands System

True precision is teamwork. The four-hands system allows both dentist and assistant to work seamlessly together, each with full freedom of movement. With this synergy, procedures flow smoothly, with efficiency and precision at every step.



Mirror Work

Vision is the heart of dentistry. **XWELLbed's** design supports both direct and reflected vision, enabling every procedure to be performed with clarity and a relaxed, natural posture - because comfort sharpens focus.




Sustainability

Perfect balance, every time


Dentistry should feel natural, never forced. That's why Swident pioneered the horizontal treatment position in Europe - a breakthrough that changed the way dentists work.

The **XWELLbed** is the evolution of this vision: a refined horizontal design that supports optimal working conditions at all times. By preserving the body's natural range of motion, it ensures fluid, fatigue free movements that protect the practitioner as much as the patient.

Left -or right- handed operation is always possible, keeping every treatment simple, flexible, and sustainable.



----- Maximum Working Area (Barnes)
-o-o- Normal Working Area (Squires)



Simple and Safe

Less is more. And safer.

With **XWELLbed**, simplicity is safety. Every detail has been carefully refined to reduce risks and extend durability.



No spittoon. Lower infection risk.

Removing the spittoon minimizes cross-contamination between patient, dentist, and assistant.



No reclining. Greater reliability.

Fewer features mean fewer failures. Every component is engineered with precision, ensuring longlasting durability. Integrated Trendelenburg movement for emergency patient positioning.



Easy maintenance. Maximum protection.

A design that's simple to clean and maintain reduces risks for staff, keeping both patients and professionals safe.

Functions

Designed for your comfort, built for your future

Every detail of **XWELLbed** is guided by ergonomic design - because a healthier posture today means a longer, more fulfilling career tomorrow.

Smart, intuitive functions make each treatment smoother and more effortless, adapting seamlessly whether you work left- or right-handed.



Headrest Memory Foam (standard)

The memory foam headrest gently embraces the patient's head, offering secure comfort while naturally aligning at the fourth cervical vertebra.



Foot Control

Full chair and instrument control at your feet - for seamless, hands-free operation.



Left- or right-handed operation

Left or right handed operation is always possible, keeping every treatment simple, flexible, and sustainable.







Headrest Gelness (optional)

The **medical gel** headrest, with its V-shape, provides comfort and stability by adapting to the neck, reducing pressure, and maintaining a neutral temperature, ensuring constant well-being and a better experience for both patient and professional.

XWELLbed

TECHNICAL DATA		
Minimum height	mm	520
Maximum height	mm	890
Width	mm	720
Length	mm	1880
Net weight approx.	Kg	130
Maximum lifting capability	Kg	200
Trendelenburg	°	20°



XWELLbed

 SW SWIDENT SA

Via G. Cattori, 3
6902 Paradiso
Lugano
Switzerland

www.swident.ch



SWIDENT SRL

Via Rossini, 26
40055 Villanova di Castenaso (BO), Italy
Tel. +39 051 6258657
info@swident.it

www.swident.it